

NewFlash Protein AnyKD PAGE Gel

Cat#: 8012011

[Product name]

NewFlash Protein AnyKD PAGE Gel

[Model & Size]

Product Name	Cat No.	Size	Component No.	Component Name	Component Size	Quantity
NewFlash Protein AnyKD PAGE Gel	8012011	50 T/Kit	8012011-1	Stacking Gel Solution A	50 mL	1 bottle
			8012011-2	Stacking Gel Solution B	50 mL	1 bottle
			8012011-3	Resolving Gel Solution A	125 mL	1 bottle
			8012011-4	Resolving Gel Solution B	125 mL	1 bottle
			8012011-5	Ammonium Persulfate (APS)	0.5 g	1 bottle

[Product Description]

NewFlash Protein AnyKD PAGE Gel kit is a product with a low acrylamide concentration, and adopts the latest gel technology to reduce the concentration of acrylamide and enhance the gelation rate. The kit is easy to use, and prepares gel quickly and safely.

Protein gel preparation is simple, and can be completed in 25 min by mixing resolving gels and stacking gels at the same time.

Protein gel electrophoresis can be completed either at 300 V constant high voltage in about 25 min or at normal voltage.

The kit is suitable for separation and identification of 10-250 kDa protein, and there is no need to adjust the concentration of resolving gel according to protein molecular weight.

[Storage and Transportation]

Stored at 2°C~8°C protected from light, with a shelf life of 24 months.

Transported on blue ice.

[The number of gels prepared per kit]

Gel Size	0.75 mm Mini-gel	1.0 mm Mini-gel	1.5 mm Mini-gel
Gel Number	62 pieces	50 pieces	33 pieces

[Operating Instructions]

1. With NewFlash Protein AnyKD PAGE Gel kit, uniformly mix resolving gel solutions A and B at a ratio of 1:1 (As a reference, for a 0.75/1.0/1.5 mm gel, respectively draw 2.0/2.5/3.8 mL of resolving gel solutions A and B). Gel preparation can be completed in a 15 mL or 50 mL centrifuge tube.

2. With NewFlash Protein AnyKD PAGE Gel kit, uniformly mix stacking gel solutions A and B at a ratio of 1:1 (As a reference, for a 0.75/1.0/1.5 mm gel, respectively draw 0.8/1.0/1.5 mL of stacking gel solutions A and B).

Note: With this product, gelation is rapid. Do not wait for pouring the resolving gel into the glass plate to start the preparation of stacking gel. Excessive waiting time will lead to partial polymerization of the resolving gel, resulting in uneven division surface of resolving gel and stacking gel, which will influence the

electrophoresis effect.

3. Add 10% APS solution to the resolving gel solution in step 1 (for 5 mL of resolving gel solution, add 50 μ L of 10% APS solution), and pour the uniformly mixed resolving gel solution into the glass plate to a position 1.5 cm from the top of the front glass plate or 0.5 cm from teeth of a comb.

Note: When it is necessary to prepare more than one protein gels at one time, the amount of APS can be appropriately reduced to lower the gelation rate. Excessive APS will result in gel embrittlement.

4. Add 10% APS solution to the stacking gel solution in step 2 (for 2 mL of stacking gel solution, add 20 μ L of 10% APS solution), and directly pour the uniformly mixed stacking gel solution on the resolving gel solution without waiting for solidification of the resolving gel solution.

Note: The kit can complete protein gel casting at one time without mixing stacking gel with resolving gel.

5. Insert the comb into the gel, and keep it standing for 15-20 min to wait for gel polymerization.

Note: Gels can be stored for several weeks at 2~8°C in sealing bags with a small amount of electrophoresis buffer.

6. After gel polymerization, remove the comb, and purge gel pores with a 1 mL injector or a pipette for sample loading.

7. The electrophoresis of this product can be carried out rapidly in a conventional electrophoresis buffer. The voltage can be set to a maximum of 300 V, at which electrophoresis can be completed in about 25 min. The current shall not exceed 140 mA. The voltage can be reduced if the current is too high (to avoid much heat).

Note: The electrophoresis of this product can be carried out at normal voltage.









[Precautions]

1. This product contains a small amount of acrylamide which is corrosive. For your safety and health, please wear a lab coat and disposable gloves.

2. This product is only used for scientific research by professionals and shall neither be used for clinical diagnosis, treatment, food or drugs nor stored in ordinary housing.

3. **Add 5 mL of deionized water to the ammonium persulfate (APS) tube to prepare 10% APS solution, and repack the solution into small portions for storage.**

[Description of Product Symbol]

Product Symbol	Description	Product Symbol	Description
	Catalog Number		Batch Code
	Date of Manufacture		Manufacturer
	Keep away from light		Temperature limit
	Consult instructions for use		Use-by date

[Company Information]

Manufacturer and after-sales service unit Name: Shenzhen Dakewe Bio-engineering Co., Ltd.

Website: www.dakewe.com

Telephone: (86-755) 86235300

Email: RD@dakewe.com

Address: Room 702-703, Building No.1, Shenzhen Biomedicine Innovations Industrial Park, No.14 Jinhui Road, Kengzi Street, Pingshan District, Shenzhen, China

After-sales service telephone: (86-755) 86235300

Zip Code: 518122