

PHA

Cat#:2030411

[Product Description]

Phytohemagglutinin isolated from *Phaseolus vulgaris* is a mixture of glycoproteins isolated from red kidney beans and contains all five isolectins found in the red kidney bean. It has been widely used for activation of lymphoid cells.

[Preparation Instructions]

PHA is supplied as a white pancake-like powder. After reconstitution is a clear solution. Do not use if the solution is or has become turbid.

Reconstitute with 500 μ L of sterile PBS before used. Stand for 5 minutes and mix gently to ensure complete solubilization. Add 10 μ L of the PHA resuspension product (25 μ g/mL) to the wells which contains 100 μ L cell culture fluid at the working concentration of 2.5 μ g/mL, which could make 10^5 PBMC produce obvious positive blots.

[Stability, Storage and Specifications]

1. Lyophilized powder is stable for 24 months from date of receipt when stored at $-25^{\circ}\text{C} \sim -15^{\circ}\text{C}$. Solution may be stored frozen in aliquots at -20°C to -70°C . PHA is one-time-use. Reuse and repeated freezing are not recommended.
2. This product is for ELISPOT quality control material or ELISPOT positive control use only, not for drug, household, or other use.
3. PHA lyophilized powder is from a 0.22 μ m filtered solution in PBS with 1%BSA.Preservative free.
4. Safety: No toxic or hazardous material.

[Company Information]

Manufacturer and after-sales service unit Name: Shenzhen Dakewe Bio-engineering Co., Ltd.

Website:www.dakewe.com

Telephone: (86-755) 86235300

Email: RD@dakewe.com

Address: Room 702-703, Building No.1, Shenzhen Biomedicine Innovations Industrial Park, No.14 Jinhui Road, Kengzi Street, Pingshan District, Shenzhen, China

After-sales service telephone: (86-755) 86235300

Zip Code: 518122