

SuperDissol Instant TBS-T Buffer

Cat#: 8015131/8015132

[Product Name]

SuperDissol Instant TBS-T Buffer

[Model & Size]

Product Name	Cat No.	Size
SuperDissol Instant TBS-T Buffer	8015131	20 packs/box
	8015132	100 packs/box

Note: The pH value of 1× solution is 7.4±0.2@25°C.

[Product Description]

Instant TBS-T Buffer is white to off-white instant granules, and each pack of this product can produce 1 L of 1 × TBS-T Buffer, which is easy to operate and use. The main components of TBS-T are sodium chloride, potassium chloride, Tris buffer salt and Tween-20. The concentration of sodium chloride, potassium chloride and Tris buffer salt in 1 × TBS-T Buffer is 137 mM, 2.7 mM and 25 mM respectively, and the content of Tween-20 is 0.05%.

Instant TBS-T Buffer is widely used in biology, which is mainly used in immunohistochemistry, in situ hybridization, ELISA and western blot, and is also commonly used as wash buffer, as well as for antibody dilution and preparation of blocking buffer.

[Storage And Transportation]

Store at room temperature of 10°C~30°C, protected from light, with a shelf life of 36 months.

[Operating Instruction]

1. Place a clean beaker, and add the entire contents of a newly opened pack of Instant TBS-T Buffer.
2. Measure about 500 ml of pure water, slowly add it into the beaker, and gently stir with a stir bar until it is completely dissolved (This step can also be performed with magnetic stirrer).
3. Transfer the TBS-T solution in Step 2 into a constant volume container, and add water to dilute to 1 L, which is 1 × working solution.

[Schematic Diagram of Instant Buffer Preparation]












Figure 1. Schematic Diagram of Instant Buffer Preparation

Note: Wash the container sufficiently with pure water before preparation, and clean it in time after preparation.

[Description Of Product Symbol]

Product Symbol	Description	Product Symbol	Description
	Catalog Number		Batch Code
	Date of Manufacture		Manufacturer
	Keep away from light		Temperature limit
	Consult instructions for use		Use-by date

[Company Information]

Manufacturer and after-sales service unit Name: Shenzhen Dakewe Bio-engineering Co., Ltd.

Website: www.dakewe.com

Telephone: (86-755) 86235300

Email: RD@dakewe.com

Address: Room 702-703, Building No.1, Shenzhen Biomedicine Innovations Industrial Park,
No.14 Jinhui Road, Kengzi Street, Pingshan District, Shenzhen, China

After-sales service telephone: (86-755) 86235300

Zip Code: 518122