

## SUPERCULTURE N-2 Supplement

Cat#: 6065011

### [Product Name]

N-2 Supplement

### [Product Description]

N-2 Supplement is a chemically determined supplement, which is mainly recommended for the growth and expression of neuroblastoma and neurons in the late stage of mitosis derived from the peripheral nervous system (PNS) and central nervous system (CNS).

This product is supplied as a 100× liquid and is intended to be used with Neuronal basal medium or DMEM/F12 medium supplemented with growth factors such as bFGF and EGF.

### [Model & Specification]

Product Name	Cat. No.	Size
N-2 Supplement	6065011	5 mL/bottle

### [Storage Conditions and Validity Period]

N-2 Supplement : Keep away from light at -15°C~-25°C, valid for 18 months.

### [Directions for Use]








Add to the medium at a final concentration of 1% (v/v).

### [Note]

1. If the supplements in this product cannot be used up at one time, it can be frozen in separate packages to avoid repeated freezing and thawing, and used within the validity period of the product. It is recommended to keep it in 2°C~8°C and keep it away from light after mixing into complete culture medium, and use it up within two weeks, and pay attention to the aseptic operation when using.

2. This supplement is a GMP-grade product for subsequent manufacturing and laboratory use.

### [Description of Product Symbol]

Product Symbol	Description	Product Symbol	Description
	Catalogue number		Batch code
	Date of manufacture		Manufacturer
	Use-by date		Temperature limit
	Consult instructions for use	/	/

### [Related Products]

Product Name	Cat. No.	Size
Neuronal Basal Medium	6016131/6016141	500 mL/bottle
NS-27 Supplement	6065021	10 mL/ bottle
ITS Supplement	6065031	10 mL/ bottle

**[Instruction Revision Date]**  
2024.03.29

**[Company Information]**

**Manufacturer and after-sales service unit Name:** Shenzhen Dakewe Bio-engineering Co., Ltd.

**Website:** www.dakewe.com

**Telephone:** (86-755) 86235300

**Email:** RD@dakewe.com

**Address:** Room 702-703, Building No.1, Shenzhen Biomedicine Innovations Industrial Park,  
No.14 Jinhui Road, Kengzi Street, Pingshan District, Shenzhen, China

**After-sales service telephone:** (86-755) 86235300

**Zip Code:** 518122