

Nylon Net

Cat#: 7061011

[Product Name]

Nylon Net

[Model & Specification]

Product Name	Cat No.	Specification	Size	Mesh Number
Nylon Net	7061011	25 Sheets/Package	90 mm×90 mm	200




[Intended Use]

Designed for use in conjunction with mouse lymphocyte separation medium. Intended for the grinding and separation of lymphocytes from the spleen of mice/rats.

[Notes]

1. This product has not undergone any sterilization process.
2. Sterilize by autoclaving at high pressure (121°C, 0.1 MPa, 30 min) before conducting aseptic procedures.
3. This product is for single use only. Please discard after use and do not reuse.

[Description of Product Symbol]

Product Symbol	Description	Product Symbol	Description
	Catalogue Number		Product Batch Code
	Consult instructions for use		

[Related Products]

Product Name	Cat No.	Size
Mouse Lymphocyte Separation Medium	7211011	100 mL/bottle
Mouse Lymphocyte Separation Medium	7211012	250 mL/bottle

[Instruction Revision Date]

2024.03.22

[Company Information]

Manufacturer and after-sales service unit Name:

Shenzhen Dakewe Bio-engineering Co., Ltd.

Website: www.dakewe.com

Telephone: (86-755) 86235300

Email: RD@dakewe.com

Production Address: Room 702-703, Building No.1, Shenzhen Biomedicine Innovations Industrial Park, No.14 Jinhui Road, Kengzi Street, Pingshan District, Shenzhen

After-sales service telephone: (86-755) 86235300

Zip Code: 518122