



Certificate of Analysis

L500 Serum-Free Medium for Lymphocyte

Cat.Number 6811021
Lot.Number A6269586
Volume 1000ml/Bottle

TEST:

TEST	REFERENCE RANGE	RESULTS
Appearance	Orange-red liquid, clear without precipitation	Pass
pH	7.10-7.50	7.41
Osmolality	300-340 mOsm/kg	328
Endotoxin	≤ 0.3 EU/mL	Pass
Sterility Testing	Negative	Pass
Mycoplasma Testing	Negative	Pass
Cell Culture Testing	Pass	Pass

RESULTS:

The product lot has met all specifications for release.



QC : Sherry

Date: Sept 14, 2023

Shenzhen Dakewe Bio-engineering Co., Ltd.

Telephone: (86-755) 86235300

Registration and Production address:

Room 702-703, Building No.1, Shenzhen Biomedicine Innovations Industrial Park, No.14
Jinhui Road, Kengzi Street, Pingshan District, Shenzhen, China

After-sales service telephone: (86-755) 86235300

Zip code: 518122