

SUPERCULTURE Cell Adhesion Reagent

Cat#: 6063011

[Product Name]

Cell Adhesion Reagent

[Product Description]

Cell Adhesion Reagent is used for tissue and cell culture to improve cell climbing efficiency in primary cell culture and cell proliferation efficiency during passaging.

[Model & Specification]

Product Name	Cat. No.	Size	Quantity
Cell Adhesion Reagent	6063011	100 mL	1 bottle

[Storage Conditions and Validity Period]

Store at -15°C~-25°C, valid for 1 year.

Store at 2°C~8°C, valid for 2 weeks.

It is recommended to aliquot and freeze the product upon arrival.








[Directions for Use]

1. Add the Cell Adhesion Reagent to the culture flask, shake the flask to evenly coat the bottom (recommended volumes: T25, 2 mL; T75, 6 mL, and so on).
2. Allow the culture flask to stand at room temperature for 1 hour or overnight at 2°C~8°C.
3. After removing the Cell Adhesion Reagent, open the culture flask cap on the clean workbench and let it stand for 15 minutes to air-dry the remaining liquid in the flask.
4. Processed culture flask can be stored at 2°C~8°C and used within a week.
5. Before cell seeding, rinse once with basic culture medium to remove residual Cell Adhesion Reagent.

[Notes]

1. This product should be used by professional personnel. Before use, the user manual should be carefully read and personal hygiene protection should be taken.
2. Pay attention to aseptic operation.
3. It is not recommended to reuse processed flask.

[Product Description of Labels and Symbols]

Product Symbol	Description	Product Symbol	Description
	Catalogue number		Batch code
	Date of manufacture		Manufacturer
	Use-by date		Temperature limitation
	Consult instructions for use	/	/

[Instruction Revision Date]

2024.03.27

[Company Information]**Manufacturer and after-sales service unit Name:**

Shenzhen Dakewe Bio-engineering Co., Ltd.

Website: www.dakewe.com

Telephone: (86-755) 86235300

Email: RD@dakewe.com

Address: Room 702-703, Building No.1, Shenzhen Biomedicine Innovations Industrial Park, No.14 Jinhui Road, Kengzi Street, Pingshan District, Shenzhen, China

After-sales service telephone: (86-755) 86235300

Zip Code: 518122