

## SUPERCULTURE NS-27 Supplement

Cat#: 6065021

### [Product Name]

NS-27 Supplement

### [Product Description]

NS-27 Supplement is the most commonly used neuron culture supplement, an optimized serum-free supplement used to support the growth and viability maintenance of embryonic, postnatal, adult hippocampal neurons and other central nervous system (CNS) neurons.

NS-27 Supplement is provided as a 50× liquid and is intended to be used with Neuronal Basal Medium for neurons culture without the need for an astrocyte feeder layers.

### [Model & Specification]

Product Name	Cat. No.	Size
NS-27 Supplement	6065021	10 mL/bottle

### [Storage Conditions and Validity Period]

NS-27 Supplement : Keep away from light at -15°C~-25°C, valid for 18 months.

### [Directions for Use]








Add to the medium at a final concentration of 2% (v/v).

### [Note]

1. If the supplements in this product cannot be used up at one time, it can be frozen in separate packages to avoid repeated freezing and thawing, and used within the validity period of the product. It is recommended to keep it in 2°C~8°C and keep it away from light after mixing into complete culture medium, and use it up within two weeks, and pay attention to the aseptic operation when using.

2. This supplement is a GMP-grade product for subsequent manufacturing and laboratory use.

### [Description of Product Symbol]

Product Symbol	Description	Product Symbol	Description
	Catalogue number		Batch code
	Date of manufacture		Manufacturer
	Use-by date		Temperature limit
	Consult instructions for use	/	/

### [Related Products]

Product Name	Cat. No.	Size
Neuronal Basal Medium	6016131/6016141	500 mL/bottle
N-2 Supplement	6065011	5 mL/ bottle
ITS Supplement	6065031	10 mL/ bottle

**[Instruction Revision Date]**

2024.03.29

**[Company Information]**

**Manufacturer and after-sales service unit Name:** Shenzhen Dakewe Bio-engineering Co., Ltd.

**Website:** www.dakewe.com

**Telephone:** (86-755) 86235300

**Email:** RD@dakewe.com

**Address:** Room 702-703, Building No.1, Shenzhen Biomedicine Innovations Industrial Park, No.14 Jinhui Road, Kengzi Street, Pingshan District, Shenzhen, China

**After-sales service telephone:** (86-755) 86235300

**Zip Code:** 518122