

## SUPERCULTURE ITS Supplement

Cat#: 6065031

### [Product Name]

ITS Supplement

### [Product Description]

ITS Supplement is a mixture of recombinant human insulin, human transferrin and sodium selenite. It is a universal supplement that can be used in conjunction with other basal media to reduce the amount of fetal bovine serum required to culture cells.

This product is provided as a 100× liquid and is intended to be used with RPMI 1640, MEM, DMEM/F12 and other basal medium.

### [Model & Specification]

Product Name	Cat. No.	Size
ITS Supplement	6065031	10 mL/bottle

### [Storage Conditions and Validity Period]

ITS Supplement : Keep away from light at -15°C~-25°C, valid for 18 months.

### [Directions for Use]








Add to the medium at a final concentration of 1% (v/v).

### [Note]

1. If the supplements in this product cannot be used up at one time, it can be frozen in separate packages to avoid repeated freezing and thawing, and used within the validity period of the product. It is recommended to keep it in 2°C~8°C and keep it away from light after mixing into complete culture medium, and use it up within two weeks, and pay attention to the aseptic operation when using.

2. This supplement is a GMP-grade product for subsequent manufacturing and laboratory use.

### [Description of Product Symbol]

Product Symbol	Description	Product Symbol	Description
	Catalogue number		Batch code
	Date of manufacture		Manufacturer
	Use-by date		Temperature limit
	Consult instructions for use	/	/

### [Related Products]

Product Name	Cat. No.	Size
RPMI-1640 Medium	6016011	500 mL/bottle
DMEM Medium -Low glucose	6016111	500 mL/ bottle
DMEM Medium -High glucose	6065121	500 mL/ bottle

**[Instruction Revision Date]**

2024.03.29

**[Company Information]****Manufacturer and after-sales service unit Name:** Shenzhen Dakewe Bio-engineering Co., Ltd.**Website:** www.dakewe.com**Telephone:** (86-755) 86235300**Email:** RD@dakewe.com**Address:** Room 702-703, Building No.1, Shenzhen Biomedicine Innovations Industrial Park,  
No.14 Jinhui Road, Kengzi Street, Pingshan District, Shenzhen, China**After-sales service telephone:** (86-755) 86235300**Zip Code:** 518122