

SUPERCULTURE DPBS

Cat#: 6062011

[Product Name]

DPBS

[Product Description]

DPBS, called “Dulbecco's Phosphate Buffered Saline”, is a balanced salt solution widely used in biology.

This DPBS does not contain calcium and magnesium, and can be used for washing cells before dissociation, diluting cells for counting, preparing reagents, etc.

[Model & Specification]

Product Name	Cat. No.	Size
DPBS	6062011	500 mL/bottle

[Specifications]

Concentration	1×
Phenol Red	Without phenol red
Calcium	Without calcium
Magnesium	Without magnesium









[Storage Conditions and Validity Period]

Stored at room temperature, the shelf life is 24 months.

[Notes]

Pay attention to operate in an aseptic condition.

[Description of Product Symbol]

Product Symbol	Description	Product Symbol	Description
	Catalogue number		Batch code
	Date of manufacture		Use-by date
	Manufacturer		Temperature limit
	Keep away from sunlight		Consult instructions for use

[Formulation]

Components	Concentration (g/L)
Potassium Chloride (KCl)	0.2
Potassium Phosphate monobasic (KH ₂ PO ₄)	0.2
Sodium Chloride (NaCl)	8
Sodium Phosphate dibasic (Na ₂ HPO ₄)	1.15

[Related Products]

Product Name	Cat. No.	Size
DMEM Medium-Low Glucose	6016111	500 mL/bottle
DMEM Medium-High Glucose	6016121	500 mL/bottle
α -MEM Medium	6016211/6016221	500 mL/bottle
IMDM Medium	6016311	500 mL/bottle
DMEM/F-12 Medium	6016021	500 mL/bottle
RPMI-1640 Medium	6016011	500 mL/bottle
Fetal Bovine Serum (South America)	6021021	500 mL/bottle
Fetal Bovine Serum (New Zealand)	6021031	500 mL/bottle
Fetal Bovine Serum (China)	6021061	500 mL/bottle
Exosome-Free Fetal Bovine Serum	6021041	50 mL/bottle

[Instruction Revision Date]

2024.03.22

[Company Information]
Manufacturer and after-sales service unit Name:

Shenzhen Dakewe Bio-engineering Co., Ltd.

Website: www.dakewe.com

Telephone: (86-755) 86235300

Email: RD@dakewe.com

Production Address: Room 702-703, Building No.1, Shenzhen Biomedicine Innovations Industrial Park, No.14 Jinhui Road, Kengzi Street, Pingshan District, Shenzhen

After-sales service telephone: (86-755) 86235300

Zip Code: 518122